



> LB Broth Miller's

<i>cat. no.</i>	<i>amount</i>	<i>note</i>
STS-LBM	1Kg	powder

For propagation and maintenance of E. Coli..Dissolve 25 gm in 1 L distilled water and adjust the pH to 7.2. Sterilize by autoclaving for 15 minutes.

Casein Digest Peptone. 10g/L

Yeast Extract 5g/L

Sodium Chloride 10g/L

Solubility

Soluble in water.

Storage

Store at room temperature.

FOR RESEARCH USE ONLY