

GeneSpin Srl Via Friuli, 51 -20135 Milano P.IVA e C.F. 04520270960 administration@pec.genespin.com info@genespin.com

>Acrylamide/ bis-Acrylamide 29.0:1 Solution

Acrylamide /bis-Acrylamide, 30% (w/v) solution

cat. no.	amount	note
STS-ACBIS 291-30	500ml	Filtered; Electrophoresis-grade.

Use this **3**0% acrylamide/bis-acrylamide, **2**9.0 :1 solution as a faster and safer alternative to handling powdered acrylamide and bis-acrylamide.

Ready-to-use high-purity (99.9%) solution

Reduce inhalation and contact hazards associated with weighing and preparing acrylamide and bis-acrylamide solutions

Precise composition and high-purity provide uniformity of gel matrices, consistent polymerization, and run-to-run reproducibility

Poison.

Acrylamide is a neurotoxin and a suspected carcinogen. Avoid contact with eyes, skin, and clothing. Wash thoroughly with soap and water after handling. See MSDS.

Store at 4°C. Keep tightly sealed. If high resolution is required (such as sequencing), fresh solution (less than 2 months) should be used.