



GeneSpin Srl Via Friuli, 51 -20135 Milano
P.IVA e C.F. 04520270960
administration@pec.genespin.com
info@genespin.com

> Protease Inhibitor Cocktail II

Protease Inhibitor Cocktail II bacterial cell

Lyophilized mixture of five protease inhibitors for use with bacterial cell extracts. Reconstitute with 1ml DMSO and 4ml water to obtain a 5ml solution. Will treat proteases extracted from 20g E.Coli/

Specificity:

Serine Proteases, Aminopeptidases,
Cysteine Proteases, Metalloproteases and
Aspartic Proteases.

Contains:

AEBSF
EDTA
Bestatin
E-64
Pepstatin A

Storage

Store at -20°C

FOR RESEARCH USE ONLY