

GeneSpin Srl Via Friuli, 51 -20135 Milano  
P.IVA e C.F. 04520270960  
administration@pec.genespin.com  
info@genespin.com

## DATA SHEET


info@genespin.com - www.genespin.com

# > bis-Acrylamide, 2% Solution

## Bisacrylamide (N,N'-Methylenebisacrylamide), 2% (w/v) solution

<i>cat. no.</i>	<i>amount</i>	<i>note</i>
STS-ACBIS 500	100ml	Filtered; Electrophoresis-grade.

### Formulation



### Molecular weight

154.17

### Quality

Filtered; Electrophoresis-grade.

Bisacrylamide is the most frequently used crosslinking agent for polyacrylamide gels.

### Solubility

Soluble in water.

### Storage

Store at 4°C. Keep tightly sealed.

FOR RESEARCH USE ONLY