



GeneSpin Srl Via Friuli, 51 -20135 Milano
P.IVA e C.F. 04520270960
administration@pec.genespin.com
info@genespin.com

> TAE 50X Solution

<i>cat. no.</i>	<i>amount</i>	<i>note</i>
STS-TAE50	5L	liquid

Solution contains:

Tris (hydroxymethyl) aminomethane: 2 M

Acetic Acid: 1 M

EDTA, Disodium Salt Dihydrate: 50 mM

Formulation

-

Molecular weight

-

Solubility

Soluble in water

Storage

Store at room temperature.

FOR RESEARCH USE ONLY