



GeneSpin Srl Via Friuli, 51 -20135 Milano
P.IVA e C.F. 04520270960
administration@pec.genespin.com
info@genespin.com

> TBS 10X Solution

<i>cat. no.</i>	<i>amount</i>	<i>note</i>
STS-TBS1	1L	liquid
STS-TBS4	4L	liquid

1X TBS solution contains 0.025M TRIS Base and 0.8% NaCl. pH (@25°C): 7.25-7.55.

Formulation

-

Molecular weight

-

Solubility

Soluble in water

Storage

Store dry at room temperature.

FOR RESEARCH USE ONLY