



genespin

DATA SHEET

info@genespin.com - www.genespin.com

## > Dextran

<i>cat. no.</i>	<i>amount</i>	<i>note</i>
STS-DEX 10	10g	powder
STS-DEX 50	50g	powder
STS-DEX 100	100g	powder

Produced from Dextran 500,000 MW. Tested for use in nucleic acid hybridizations.  
Purity :>98%

Formulation

-

Molecular weight

-

Solubility

Soluble in water.

Storage

Store dry at room temperature.

FOR RESEARCH USE ONLY