



genespin

DATA SHEET

info@genespin.com - www.genespin.com

> Nutrient Agar

<i>cat. no.</i>	<i>amount</i>	<i>note</i>
STS-NAG 500	500g	powder
STS-NAG 1000	1Kg	powder

For cultivating a variety of non-fastidious organisms. Dissolve 23g in one liter distilled water. Heat solution to boiling with continuous stirring for one minute. Sterilize by autoclaving at 121°C for 15 minutes. Mix medium with gentle agitation to avoid bubbles prior to pouring plates.

Beef Extract 3 g/l
Pancreatic Digest of Gelatin 5 g/l
Agar 15 g/l

Formulation

-

Molecular weight

-

Solubility

Soluble in water.

Storage

Store at room temperature.

FOR RESEARCH USE ONLY