



genespin

DATA SHEET

info@genespin.com - www.genespin.com

protease

> PBS 10X Solution

<i>cat. no.</i>	<i>amount</i>	<i>note</i>
STS-PBS1	1L	liquid
STS-PBS4	4L	liquid

PBS is a common buffer used in many laboratory procedures. PBS helps to maintain a constant pH.

Solution contains:

Sodium Chloride 8.0%

Potassium Chloride <1.0%

Sodium Phosphate Dibasic 1.4%

Potassium Phosphate Monobasic <1.0%

Formulation

-

Molecular weight

-

Solubility

Soluble in water

Storage

Store dry at room temperature.

FOR RESEARCH USE ONLY