



genespin

DATA SHEET

info@genespin.com - www.genespin.com

> TBE 5X Solution

<i>cat. no.</i>	<i>amount</i>	<i>note</i>
STS-TBE 1	1L	5x conc. solution
STS-TBE 4	4L	5x conc. solution

Formulation

Dry homogenous blend of molecular biology grade Tris-base boric acid and EDTA, disodium salt. For use in gel electrophoresis.

Solubility

Soluble in water.

Storage

Store at room temperature.

FOR RESEARCH USE ONLY