

GeneSpin Srl Via Friuli, 51 -20135 Milano P.IVA e C.F. 04520270960 administration@pec.genespin.com info@genespin.com

## > PBS 100mL TABLETS

cat. no.	amount	note
STS-PBST	100	tablets

Phosphate Buffered Saline (PBS) Tablets. Add one tablet to 100ml of purified water for a 1X solution of PBS. Contains 137mM Sodium Chloride, 2.7mM Potassium Chloride and 10.0mM Phosphate.

Formulation

Molecular weight

Solubility Soluble in water

Storage Store dry at room temperature.