



GeneSpin Srl Via Friuli, 51 -20135 Milano  
P.IVA e C.F. 04520270960  
administration@pec.genespin.com  
info@genespin.com

**phosphatase**

## > Phosphatase Inhibitor Cocktail II

Phosphatase Inhibitor Cocktail II

Lyophilized blend of five phosphatase inhibitors for use with animal tissue. A431 and Jurkat cell extracts. Suitable for extracts containing detergents. Add 1ml of pure water to produce a 100x concentrated. Dilute 1:100 before using.

Specificity:

Tyrosine Phosphatases, Acid Phosphatases and Alkaline Phosphatases

Contains:

Imidazole, Sodium Fluoride,  
Sodium Molybdate, Sodium  
Orthovanadate, Sodium Tartrate,  
dihydrate

Storage

Store at -20°C

FOR RESEARCH USE ONLY