

GeneSpin Srl Via Friuli, 51 -20135 Milano P.IVA e C.F. 04520270960 administration@pec.genespin.com info@genespin.com

## > Phosphatase Inhibitor Cocktail III

Phosphatase Inhibitor Cocktail III

Lyophilized blend of four phosphatase inhibitors for use with animal tissue estracts. Reconstitute with 1ml water to produce a 1 ml concentrate. Dilute 1:5 before using.

Specificity: Ser/Thr Phosphatases and Protein Tyrosine Phosphatases

Contains:
Sodium Fluoride,
Sodium Orthovanadate
Sodium Phyrophosphate, decahydrate
B-Glycerophosphate

Storage Store at -20°C