



**phosphatase**

## > Phosphatase Inhibitor Cocktail III

Phosphatase Inhibitor Cocktail III

Lyophilized blend of four phosphatase inhibitors for use with animal tissue extracts. Reconstitute with 1ml water to produce a 1 ml concentrate. Dilute 1:5 before using.

Specificity:

Ser/Thr Phosphatases and Protein Tyrosine  
Phosphatases

Contains:

Sodium Fluoride,  
Sodium Orthovanadate  
Sodium Phyrophosphate, decahydrate  
B-Glycerophosphate

Storage

Store at -20°C

FOR RESEARCH USE ONLY