



GeneSpin Srl Via Friuli, 51 -20135 Milano  
P.IVA e C.F. 04520270960  
administration@pec.genespin.com  
info@genespin.com

## > Phosphoramidon

cat. no.	amount	note
STS-PHR 5	5mg	powder

Purity: >90%

Protease inhibitor specific for thermolysin and other metalloproteases. Also effective inhibitor of collagenase. Working Concentration: 5-500um

Formulation

$C_{23}H_{34}N_3O_{10}P$

Molecular weight

543.5

Solubility

Soluble in DMSO and methanol

Storage

Store dry at -20°C.

FOR RESEARCH USE ONLY