



DATA SHEET

GeneSpin Srl Via Friuli, 51 -20135 Milano  
P.IVA e C.F. 04520270960  
administration@pec.genespin.com  
info@genespin.com

**protease**

## > Protease Inhibitor Cocktail I

Protease Inhibitor Cocktail I for general use

Lyophilized blend of five inhibitors for inhibiting a broad range of proteases and esterases. Reconstitutes with 1 ml of water to produce a 100X concentrate solution.

Specificity:

Serine Proteases, Esterases, Cysteine  
Proteases, Metalloproteases  
and Trypsin like Protease

Contains:

AEBSF E-64 Leupeptin  
Aprotinin EDTA

Storage

Store at -20°C

FOR RESEARCH USE ONLY