



DATA SHEET

protease

> Protease Inhibitor Cocktail I

Protease Inhibitor Cocktail I for general use

Lyophilized blend of five inhibitors for inhibiting a broad range of proteases and esterases. Reconstitutes with 1 ml of water to produce a 100X concentrate solution.

Specificity:

Serine Proteases, Esterases, Cysteine
Proteases, Metalloproteases
and Trypsin like Protease

Contains:

AEBSF E-64 Leupeptin
Aprotinin EDTA

Storage

Store at -20°C

FOR RESEARCH USE ONLY