



GeneSpin Srl Via Friuli, 51 -20135 Milano
P.IVA e C.F. 04520270960
administration@pec.genespin.com
info@genespin.com

> Protease Inhibitor Cocktail V

Protease Inhibitor Cocktail V EDTA free

Lyophilized blend of four protease inhibitors for inhibiting serine and cysteine but not metallo proteases. Reconstitute with 1ml of water to produce a 100X stock solution

Specificity:

Serine Proteases, Cysteine Proteases, and
Trypsin like Proteases

Contains:

AEBSF E-64

Aprotinin Leupeptin

Storage

Store at -20°C

FOR RESEARCH USE ONLY