



## > Protease Inhibitor Cocktail VII

Protease Inhibitor Cocktail VII for HIS-Tag sequence

Mixture of five protease inhibitors for purification of proteins containing hexa histidine tagged sequences. Supplied in 1ml DMSO. Will inhibit proteases extracted from 10g of cells.

Specificity:

Serine Proteases, Cysteine Proteases,  
Aspartic Proteases, Aminopeptidases  
and Metalloproteases

Contains:

AEBSF  
Bestatin  
E-64  
Pepstatin A  
Phosphoramidon

Storage

Store at -20°C

FOR RESEARCH USE ONLY